Sable Industries Inc.
100 Campbell Avenue, Suite 5,
Kitchener, Ontario N2H 4X8 Canada
www.sableindustriesinc.com

SABLE SLOW SPEED E-TYPE AIR MOTOR MAINTENANCE AND CLEANING MANUAL

Please retain this manual as your proof of purchase.

MODEL NO.

- ☐ 1600101
☐ 1600102

MAXIMUM RPM

- 25,000
5,000

FEATURES

1. Air Motor conforms to International Standard Organization (ISO) specifications.
2. Accepts all E Type (ISO 3964 compatible) attachments.
3. Autoclavable.

SPECIFICATIONS

- | | |
|--------------------------------|-----------------------|
| 1. Maximum Air Pressure | : 45 PSI |
| 2. Connection with Attachments | : ISO Standard |
| 3. Method of Sterilization | : Autoclave |
| 4. Air hose connection | : Midwest type 4 hole |

HOW TO USE

Use only Clean, Dry, Compressed Air.

1. Connect to air hose. (Ensure rubber gasket is in place.)
2. Select Forward or Reverse Rotation

Caution: For maximum speed, selection ring must be fully in F or R position.

3. Attachment of E Type Contra Angles or Nosecones.

Slide Attachment onto shaft of Air Motor until locked in place.

4. Removal of Contra Angles or Nosecones.

Pull Attachment, pushing release button located on the exterior of Air Motor.

Note: Exercise caution when operating. For use only by trained personnel.

LUBRICATION

USE ONLY AUTOCLAVABLE HANDPIECE LUBRICANT

1. Remove Handpiece from air hose.
2. Apply lubricant into drive air tube Fig. 1.
3. Attach Handpiece to air hose.
4. Run Handpiece at normal speed for 10 seconds.
5. Wipe off excess oil from Handpiece surface.

AUTOCLAVING PROCEDURE

1. Wipe off debris on the Handpiece surface with alcohol-soaked cloth.
2. Lubricate the Handpiece. (Refer to Lubrication Instructions).
3. Autoclave for 20 minutes at 135 C.

ATTENTION

This is a precision instrument. Do not drop or apply excessive force when handling.

WARRANTY INFORMATION

The Sable slow speed Handpiece is warranted to be free of manufacturing defects for a period of 12 months* from date of purchase.

*Warranty is conditional on proper use and maintenance of Handpiece.

Warranty is limited to repair or replacement at manufacturers discretion. Under no circumstances will the manufacturer be held responsible for incremental or consequential damages.

Please contact point of purchase authorized distributor for warranty repair.

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